*...going one step further*

**1012852**  
(W19925)

## Scope of Supply

The model 'Mini' contains the following items:

- Electrophoresis chamber with corrosion resistant platinum electrodes and high-quality connections for plus and minus pole
- Safety lid with attached power cords
- 2 combs, 1.5mm thick, 12 teeth each
- Instruction manual

## Safety Instructions

- Please read the instruction manual carefully before using the electrophoresis system.
- Use only CE marked power supplies.
- Always disconnect the electrophoresis chamber from the power supply before removing the safety lid.
- Do not put the safety lid down on wet surfaces.
- Always disconnect the electrophoresis chamber from the power supply when not in use.
- The chamber is provided for a max. current and voltage. Do not exceed 150 mA, 120 V.
- Please mind the maximum filling level for the electrophoresis buffer.

## Technical Properties

The horizontal electrophoresis system 'Mini' has been designed in such a way that allows the user to cast and run the gel directly in a single chamber. No additional casting equipment, such as grease or seal is required to prepare the gel tray for casting the gel.

The chamber's bottom is UV-permeable, thus the gel run can easily be documented in case of using DNA-fluorescent dyes. Self evidently, the gel can also be removed from the chamber for staining and documentation after the gel run.

**Please use only UV light with a wave length of 300nm or higher for all applications, as hard UV light (<300 nm) damages the gel chamber and the nucleic acids.**

If you need only short runs are needed, you can use both combs which then gives twice as many slots for samples.

## Warranty

The manufacturer guarantees that the electrophoresis chamber was completely tested before delivery and that it complies with the applicable safety regulations.

Please check the delivery immediately upon receipt for completeness and possible damage in transit. If anything is faulty or damaged, please contact 3B Scientific immediately.

The manufacturer gives a 24-month warranty on the product conditional on the chamber being used according to the instruction manual. Claims for replacement or repair are not valid if physical abuse is the reason for the damage. Any liability for consequential damages resulting from the use of the electrophoresis chamber is hereby excluded. The manufacturer's liability for intent and gross negligence or for damages resulting from injury to life, body or health is unaffected by that.

Our experience shows the two platinum electrodes last for many years if the product is used properly. Any Damage occurring is usually due to mechanical errors (e. g. by washing-up brushes) or contact with so-called "platinum poisons" such as phosphorus, boron and heavy metals (lead, zinc, etc.). For this reason, the platinum electrodes are excluded from the manufacturer's warranty.

In order to introduce developments promptly, the manufacturer reserves the right to change specifications and small visual details without prior notice.

## General Instructions

### Casting of agarose gels

1. Please wear protective clothing, goggles and gloves for heating the agarose solution.
2. Remove the safety lid of the electrophoresis chamber carefully.
3. For the preparation of the gel solution, only suitable agaroses and electrophoresis buffers should be used. For preparing the agarose solution, weigh a suitable amount of agarose in a small Erlenmeyer flask, add an appropriate amount of electrophoresis buffer and a stir bar (if heating with heating stirrers). Note the weight of the filled flask so that any loss during heating can be compensated for by the addition of distilled H<sub>2</sub>O. Thus, the gel solution will have the agarose concentration

that you need. To dissolve the agarose, heat the Erlenmeyer flask either in a microwave or with a heating stirrer. For the latter, use medium heating power and stir constantly. Heating in a microwave should only be for short periods of time, on medium heating level, repeated multiple times. When heating in a microwave, remove the Erlenmeyer flask from time to time (remember to wear gloves, goggles and be careful to avoid boiling!) and gently swirl the gel solution in a circle. Put it back in the microwave afterwards and repeat the entire process 3-4 times, until the agarose is completely dissolved.

4. Before the casting of the gel, allow the gel solution to cool down to 60°C. Cast the gel solution without any air bubbles in the gap between the two white plastic stripes on the chamber's ground (its capacity is about 40 ml). Then insert the comb (or both combs respectively). Make sure that there are no air bubbles on the edges of the slots. Gels made of standard agarose solidify within 20 minutes at room temperature.
5. After the solidification of the agarose gel, add some of the electrophoresis buffer in 3-mm layers in the chamber over the gel to prevent it from drying out. After carefully removing the comb(s), the gel can now be loaded. The coated gel can be stored for a few days at 4°C. For best results leave the comb inserted and cover the gel with wrapping film to protect it from drying out.

### **Loading of gels and electrophoresis**

1. After the solidification of the agarose gel, cover the gel with the electrophoresis buffer. Please note the maximum filling level for the electrophoresis buffer.
2. Remove the comb(s) from the agarose gel by gently moving back and forth and then carefully pulling out upright.
3. In the next step, the slots in the gel can be loaded with the samples. Preparing the samples with an appropriate gel loading buffer which increases the specific weight of the samples makes it easier to pour them into the slots with a micropipette (see section "gel loading buffer").
4. When loading the gel slots, dip the tip of the micropipette carefully into the slot and then slowly press out the sample. Make sure not to damage the bottom of the slot.  
Beginners should practise this with „practice samples“ consisting of H<sub>2</sub>O and gel loading buffer only.
5. At least one size standard should be applied on each gel so that the size of the various DNA fragments can be determined.
6. Now place the safety lid on the electrophoresis chamber and connect the chamber to a suitable power supply. Please note the correct polarity. In an alkaline up to a neutral medium, nucleic acids are charged negatively and migrate to the anode (red pole).
7. Turn on the power supply and run the electrophoresis at a voltage in the range 80 - 120 V. The process of the electrophoresis can be observed with the help of the dye, which is inside the gel loading buffer.
8. Bromophenol blue and xylene cyanol are frequently used dyes for the gel loading buffer. These dyes are also negatively charged and migrate to the anode. This process, the so-called co-migration with the double-stranded DNA fragments, depends on different factors, such as agarose type, electrophoresis buffer and gel strength.  
For a rough estimation: bromophenol blue migrates in 1xTAE electrophoresis buffer and 1%-standard agarose gel in the same way as a DNA fragment of 650 base pairs. Under the same conditions, xylene cyanol migrates like a fragment with 5 000 base pairs.

### **Staining of nucleic acids**

Since there are different methods of staining nucleic acids, we refer at this point to the corresponding technical literature (Sambrook et. al., 1989).

### **Evaluation / documentation**

Size determinations of nucleic acid fragments can be carried out by comparison with the fragments of a length standard. A detailed description of this procedure can be found, for example, in Molecular Genetics by R. Knippers (2008).

## Cleaning and Maintenance

**Attention:** Disconnect the electrophoresis chamber from power supply before cleaning.

The gel chamber and the combs should be cleaned with warm water immediately after each use. If necessary, a few drops of a washing-up liquid can be added. Please do not use dishwashing brushes or the like, as this can damage the platinum electrodes. Rinse the electrophoresis chamber after cleaning with distilled or demineralised water to prevent the formation of limescale.

Never leave the chamber uncleaned for hours or even overnight as residues from gel or buffer can cause stains which are difficult to remove completely.

**Attention:** Please avoid the electrical connections being in contact with water. If however this does occur, please carefully wipe dry with a soft cloth and allow to air-dry. Please do not use paper because it is often too rough and causes small scratches. Do not use a hair dryer either, as the hot air can damage connectors, fittings and the acrylic material.

**Attention:** Do not use ethanol nor other organic solvents for cleaning. It can cause cracks, scratches or other unwanted material changes (blindness etc.).

## Required Material and Processes

### Electrophoresis buffer

The electrophoresis buffer provides the necessary ions for the electrophoresis process. Its addition provides a constant pH value so that nucleic acids have the desired net charge. Nucleic acids are negatively charged in an alkaline up to neutral medium. Normally, electrophoresis buffer contains components that protect nucleic acids from degradation, e.g. EDTA, which complexes divalent cations and therefore inhibits DNases.

At this point, the frequently used TAE electrophoresis buffer for electrophoresis of DNA is described for non-denaturing conditions. TAE stands for Tris-acetate-EDTA. You can either produce the buffer yourself or order it from 3B Scientific as a buffer concentrate

### Agarose, gel volumes and gel concentrations

Even though various agaroses are offered, the so called standard agarose is certainly of prime importance. For the casting of an agarose gel in this chamber about 40 ml of gel solution is needed. Please prepare a little extra gel solution to allow for the fact that small residues always remain in the vessel.

Optimal separation of molecules of different sizes will be achieved by varying the concentration of the gel according to the information in the chart below.

### Optimal separation ranges of the following DNA length (double-stranded DNA) corresponding to various concentrations of agarose (standard agarose)

concentration of agarose (%)	agarose (g)	buffer (ml)	optimal separation range (kbp)
0,5	0,25	50	1 – 15
0,7	0,35	50	0,8 – 10
1,0	0,5	50	0,5 – 7
1,2	0,6	50	0,3 – 6
1,5	0,75	50	0,2 – 4
2,0	1,0	50	0,1 – 3

### Gel loading buffer

The samples which are to be analyzed need to be mixed with a suitable gel loading buffer before applying them to the gel. Gel loading buffer contains dye(s) for tracking the run of the electrophoresis as well as glycerin, saccharose and the like. Therefore, the prepared samples are heavier than the electrophoresis buffer and will decrease slightly into the gel slots during the sample application. Bromophenol blue and xylene cyanol are frequently used dyes for gel loading buffers. The patterns of their run (or the so called co-migration to double-stranded DNA fragments) depend on the type of agarose, gel strength and the type of electrophoresis buffer.

For a rough estimation: bromophenol blue migrates in 1xTAE electrophoresis buffer and 1%-standard agarose gel in the same way as a DNA fragment of 650 base pairs. Under the same conditions, xylene cyanol migrates like a fragment with 5 000 base pairs.

### **DNA length standard**

For the sizing of the samples, apply one DNA size standard on each gel at least to one lane. DNA size markers are made of DNA fragments of familiar sizes.

### **Staining of nucleic acids**

Since there are different methods of staining nucleic acids, we refer at this point to the corresponding technical literature (Sambrook et. al., 1989).

### **Declaration of Conformity**


The electrophoresis chamber 'Mini' complies with the applicable security requirements according to EN: 50081-1 (1992) and EN: 50082-1 (1992).

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